

*Newhëëk*

---

*X Ray Exposure*

*Hand Switch*

*Professional  
X Ray Machine  
Accessories Producer*

*Website : [www.xrayswitch.com](http://www.xrayswitch.com)*

*Tel: 86 536 8882360*

*Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)*



**X-ray hand switch L01**

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)



**X-ray hand switch L02**

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)



**X-ray hand switch L03**

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)



**X-ray hand switch L04**

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)





**X-ray hand switch L05**


***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)



**X-ray hand switch L06**

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

## Quick Details

Properties: Medical X-ray Equipments & Accessories

Brand Name: Newheek

Product Name: X-ray exposure hand switch

Place of Origin: Shandong, China (Mainland)

Application: x ray machine

Housing materials: ABS plastic

Standard core: Type L01/L02/L04    3 core or 4 core    can be choosed  
Type L03                            4 core or 6 core    can be choosed  
Type L05/L06                       2 core

Standard Length: 3m wire ( 30cm spring coil wire )  
4m wire ( 40cm spring coil wire )  
6m wire ( 60cm spring coil wire )  
10m wire ( 100cm spring coil wire )  
Other requirement .....

Button: Type L01/L02/L04    two-button design  
Type L03                        three-button design  
Type L05/L06                   one-button design

Structure: omron inner switch

Customization: Available

## Product Application

### **Applicable to x-ray exposure of radiography or fluoroscopy equipment**

X ray hand exposure hand switch is mainly used on portable x ray, mobile x ray, stationary x ray, analog x ray, digital x ray, radiography x ray etc x ray equipment.

Especially it can match with Siemens , Philips or GE x ray equipment.

As well as it is applicable to beauty laser device, healthy recovery device etc field.

***Weifang Newheek Medical Image Equipment Co., Ltd.***

*Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.*

*Website: [www.xrayswitch.com](http://www.xrayswitch.com)*

*Tel: 86 536 8882360*

*Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)*



## HighLight

- Streamline design outlook, comfortable handle.
- Housing adopt ABC engineering plastic, anti-shock and fire retardant.
- Inner adopt OMRON micro switch, high sensitivity.

## Technical Features

### 3-core switch      Two-button

Housing Material	Hand Switch Cable Core		
	Red	White	Green
white, ABS plastic	I Step	Common Point	II Step
Ambient Temperature	Relative Humidity		Atmospheric Pressure
-20°C-70°C	≤ 93%		50-106kPa

### 4-core switch      Two-button

Housing Material	Hand Switch Cable Core			
	White	Black	Red	Green
white, ABS plastic	I Step		II Step	
Ambient Temperature	Relative Humidity		Atmospheric Pressure	
-20°C-70°C	≤ 93%		50-106kPa	

**Weifang Newheek Medical Image Equipment Co., Ltd.**

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

**4-core switch    Three-button**

Housing Material	Hand Switch Cable Core			
	Red	White	black	Green
white, ABS plastic	I Step	Common Point	II Step	III Step
Ambient Temperature	Relative Humidity		Atmospheric Pressure	
-20°C-70°C	≤ 93%		50-106kPa	

**6-core switch    Three-button**

Housing Material	Hand Switch Cable Core					
	Red	Pink	Blue	Purple	Black	Brown
white, ABS plastic	I Step		II Step		III Step	
Ambient Temperature	Relative Humidity			Atmospheric Pressure		
-20°C-70°C	≤ 93%			50-106kPa		

**2-core switch    One-button**

Housing Material	Hand Switch Cable Core		LED	
	Red	Green	L05	L06
white, ABS plastic	I Step		yes	no
Ambient Temperature	Relative Humidity		Atmospheric Pressure	
-20°C-70°C	≤ 93%		50-106kPa	

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

**Configuration:**

Standard	3-core, 3m-wire	4-core, 3m-wire 4-core, 6m-wire	2-core, 3m-wire	6-core, 6m-wire
Customized	3-core, 6m or 10m wire	4-core, 4m or 10m wire	2-core, 6m or 10m wire	6-core, 4m or 10m wire
	Other requirement	Other requirement	Other requirement	Other requirement

**Our Service**

**Original manufacturer of x- ray machine accessories hand switch and foot switch for more than 16 years.**

√ Customers could find all kinds of x-ray machine parts here.

√ Offer on line technological support.

√ Promise super product quality with best price and service.

√ Support the third part inspection before delivery.

√ Ensure the shortest delivery time.

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

**Contact us**

Email : [info@xrayswitch.com](mailto:info@xrayswitch.com)

Hot Line: +86 536 8882360

**FAQ**

The common problems on the **x ray exposure hand switch**

1. Can you accept customize hand switch?

Yes. Please feel free to send us your detail requirement by email here.

You will get our reply in 24 hours.

2. Can we print our logo/website/company name on products?

Yes, please advise the size and Pantone Code of the logo.

3. What's the lead time for regular order?

3-5 days for order quantity lower than 100pc hand switch; For more quantity, usually 15 days.

4. Can I get a discount?

Yes, for the order quantity more than 50 pcs hand switch, please contact us to get the best price.

5. What is the quality warranty?

We offer 12 months warranty after the hand switch is delivered.

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

Certificate EMC

Certificate of Conformity

Certificate Number: B-E16029006  
EMC Directive 2014/30/EU



Holder.....: Weifang Newheek Electronic Technology Co.,Ltd.  
Address.....: The North 2nd Floor of No.1 Factory Building of No.13426 Yuqing East Street, Hi-Tech Industrial Development Zone, Weifang, Shandong, China  
Manufacturer.....: Same As Holder  
Product.....: Hand Switch  
Model No.....: XC-01, XC-02, XC-03, XC-04, L01, L02, L03, L04, L05, L06, L07, L08, L09, L10  
Technical Data.....: AC 110-240V~, 10A / DC 30V, 5A

The submitted products have been tested by us with the following standard(s) and found to be in compliance with the listed European Directives.

- EN 61000-6-3:2007+A1:2011;
- EN 61000-6-1:2007;
- EN 61000-3-2:2006+A1:2009+A2:2009;
- EN 61000-3-3:2013

The test results apply only to the particular sample tested and to the specific tests carried out. Technical Report and documentation are at the Holder's disposal.

This certificate applies specifically to the sample investigated in our test reference number only. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation. Other relevant Directives have to be observed.



*[Signature]*

Certification Manager  
Date: Feb. 24, 2016



Company No.07113834



Beide (UK) Product Service Limited

U.K.: Flat 107, 25 Indecon Square, London, United Kingdom  
China: 6F, Bldg E, Hourui 3rd Ind Zone, Xixiang, Bao'an Dist, Shenzhen, China  
Http://www.szbeide.com E-mail: admin@szbeide.com

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.  
Website: [www.xrayswitch.com](http://www.xrayswitch.com) Tel: 86 536 8882360 Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

Certificate ROHS

Certificate of Conformity

Certificate Number: B-R16029007  
RoHS Directive 2011/65/EU



Holder.....: Weifang Newheek Electronic Technology Co.,Ltd.

Address.....: The North 2nd Floor of No.1 Factory Building of No.13426 Yuqing East Street, Hi-Tech Industrial Development Zone, Weifang, Shandong, China

Manufacturer.....: Same As Holder

Product.....: Hand Switch

Model No.....: XC-01, XC-02, XC-03, XC-04, L01, L02, L03, L04, L05, L06, L07, L08, L09, L10

The submitted products have been tested by us and found to be in compliance with the listed European Directives.

The test results apply only to the particular sample tested and to the specific tests carried out. Technical Report and documentation are at the Holder's disposal.

This certificate applies specifically to the sample investigated in our test reference number only. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation. Other relevant Directives have to be observed.



  
 Certification Manager  
 Date: Feb. 24, 2016



Company No. 07113834



Beide (UK) Product Service Limited

U.K.: Flat 107, 25 Indescon Square, London, United Kingdom  
 China: 6F, Bldg E, Hourui 3rd Ind Zone, Xixiang, Bao'an Dist, Shenzhen, China  
 Http://www.szbeide.com E-mail: admin@szbeide.com

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.  
 Website: [www.xrayswitch.com](http://www.xrayswitch.com) Tel: 86 536 8882360 Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)